

Please note: The concentration of the pUC19 DNA included with this product for use as a transformation positive control has been changed from 1ng/μl to 10pg/μl to eliminate the need for a dilution step before transformation. Do not dilute the pUC19 plasmid before performing the transformation positive control. Please contact Lucigen if you have any questions.

***E. coli*[®] EXPRESS**

Chemically Competent Cells

IMPORTANT!
-80°C Storage Required
Immediately Upon Receipt

Lucigen[®] Corporation
Advanced Products for Molecular Biology

2120 W. Greenview Drive
Middleton, WI 53562

Toll Free (888) 575-9695

Phone (608) 831-9011

FAX (608) 831-9012

www.lucigen.com

E. coli[®] EXPRESS Chemically Competent Cells

Contents

Components & Storage Conditions	3
<i>E. coli</i> EXPRESS Chemically Competent Cells.....	4
Preparation for Transformation.....	4
Transformation Protocol	4
Media Recipes.....	5
Related Products.....	5

Notice of Limited Label License, Copyright, Patents, Warranties, Disclaimers and Trademarks

Copyright© 2010 by Lucigen Corp. Lucigen, BigEasy, ClonePlex, CloneSmart, CopyRight, DNATerminator, *E. coli*, EconoTaq, Lucigen, NanoClone, PCR Terminator, pJAZZ, pSMART, pTrueBlue, and PyroPhage are registered trademarks of Lucigen Corporation. bSMART, Blue/White cloning without the blues, ChimeraFree, CloneDirect, cSMART, eLucidations, Espresso, GapFree, HI-Control, PCR SMART, pEZ, pEZSeq, pGC, pLEXX-AK, pRANGER, PyroScript, Replicator, The Cloning Company, The Molecular Cloning Company, Transformance, TSA, and UltraClone are trademarks of Lucigen Corporation. Hydroshear and GeneMachines are trademarks of Genomic Solutions, Inc. OverExpress is a trademark of Imaxio S.A. (Saint-Beauzire, France).

Lucigen's products are sold for research use only and are not to be used in humans or for medical diagnostics. Lucigen's liability with respect to any CloneSmart product is limited to the replacement of the product. No other warranties of any kind, expressed or implied, including without limitation, any implied fitness for any particular use, are provided by Lucigen. Lucigen is not liable for any direct, indirect, incidental or consequential damages arising out of or in connection with the use or inability to use any of its CloneSmart products.

Limited Label License

E. coli[®] EXPRESS products are sold under a license issued to Lucigen Corporation by Brookhaven Science Associates, LLC for noncommercial research purposes only. A separate license is required for any commercial use, including the use of these materials for research purposes or production purposes by any commercial entity. Information about commercial licenses may be obtained from the Office of Technology Commercialization and Partnership, Brookhaven National Laboratory, Bldg. 490-C, P. O. Box 5000, Upton, New York 11973-5000, telephone (631)-344-7134. In using this product, you agree that *E. coli* EXPRESS cells or its derivatives containing the cloned gene for T7 RNA Polymerase may not be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of this limited label license and agrees to be bound by its terms.

If the purchaser is not willing to accept these use limitations, Lucigen Corporation is willing to accept return of the product for a full refund. For information on obtaining a license, contact Lucigen Corporation, 2120 W. Greenview Dr., Middleton, WI 53562. Email: Lucigen@lucigen.com. Phone: 608-831-9011. Fax 608-831-9012.

E. coli[®] EXPRESS Chemically Competent Cells

Components & Storage Conditions

Three strains of Lucigen's *E. coli* EXPRESS Chemically Competent Cells are available: BL21(DE3), BL21(DE3)pLysS, and BL21(DE3)pLysE.

All three strains of *E. coli* EXPRESS Chemically Competent Cells have a transformation efficiency yield of $\geq 1 \times 10^7$ cfu/ μ g pUC19.

The cells are shipped on dry ice in one container, along with supercoiled control pUC19 DNA at 10 pg/ μ l and Expression Recovery Medium (lactose minus). *E. coli* EXPRESS Chemically Competent Cells are available in 80 μ l aliquots (DUOs), sufficient for 2 transformations per tube. Please refer to the table below for catalog numbers.

E. coli[®] EXPRESS Chemically Competent Cells: Store at -80° C

STRAIN	Efficiency (cfu/ μ g pUC19)	Transformations	Catalog #	Storage
<i>E. coli</i> EXPRESS BL21(DE3) DUOs (Orange tubes)	$\geq 1 \times 10^7$	12 (6 x 80 μ l) 24 (12 x 80 μ l) 48 (24 x 80 μ l)	60401-1 60401-2 60401-3	-80° C
<i>E. coli</i> EXPRESS BL21(DE3)pLysS DUOs (Pink tubes)	$\geq 1 \times 10^7$	12 (6 x 80 μ l) 24 (12 x 80 μ l) 48 (24 x 80 μ l)	60413-1 60413-2 60413-3	-80° C
<i>E. coli</i> EXPRESS BL21(DE3)pLysE DUOs (Purple tubes)	$\geq 1 \times 10^7$	12 (6 x 80 μ l) 24 (12 x 80 μ l) 48 (24 x 80 μ l)	60425-1 60425-2 60425-3	-80° C
<i>E. coli</i> EXPRESS BL21(DE3) ChemComboPack (4 reactions of the above 3 strains)		12 (6 x 80 μ l)	60430-1	-80° C
Expression Recovery Medium (lactose minus)		12 (1 x 12 mls) 24 (2 x 12 mls) 48 (4 x 12 mls)	----	-20 to -80° C
Supercoiled pUC19 DNA (10 pg/ μ l)		(1 x 20 μ l)	----	-20 to -80° C

* Additional Expression Recovery Medium (lactose minus) can be ordered separately as Catalog # 80030-1 (8 x 12 ml)

E. coli[®] EXPRESS Chemically Competent Cells

E. coli EXPRESS BL21(DE3), BL21(DE3)pLysS, and BL21(DE3)pLysE Chemically Competent Cells are *E. coli* strains that are ideal for routine protein expression applications. NOTE: for expressing toxic proteins, we recommend Lucigen's OverExpress™ C41(DE3) and C43(DE3) Competent Cells.

BL21(DE3), BL21(DE3)pLysS, and BL21(DE3)pLysE are lysogens of λ DE3. These strains carry a chromosomal copy of the T7 RNA Polymerase gene under the control of the *lacUV5* promoter. These strains are suitable for production of protein from target genes cloned into T7 driven expression vectors. *E. coli* Express BL21(DE3), BL21(DE3)pLysS, and BL21(DE3)pLysE are also deficient in the *lon* and *ompT* proteases.

E. coli EXPRESS BL21(DE3)pLysS, and BL21(DE3)pLysE have the same chromosomal genotype as BL21(DE3). These strains also carry a chloramphenicol resistant plasmid that encode T7 lysozyme which is a natural inhibitor of T7 RNA Polymerase. Cells containing pLysS produce a small amount of T7 Lysozyme, whereas the pLysE strain produces considerably more T7 Lysozyme. These strains are used to suppress basal expression of T7 RNA polymerase prior to induction, thus stabilizing recombinants encoding potentially deleterious proteins.

Genotypes

E. coli EXPRESS BL21(DE3)

F^- *ompT hsdS_B (r_B⁻ m_B⁻) gal dcm* (DE3)

E. coli EXPRESS BL21(DE3)pLysS

F^- *ompT hsdS_B (r_B⁻ m_B⁻) gal dcm* (DE3) pLysS (Cm^R)

E. coli EXPRESS BL21(DE3)pLysE

F^- *ompT hsdS_B (r_B⁻ m_B⁻) gal dcm* (DE3) pLysE (Cm^R)

As a control for transformation, *E. coli* EXPRESS Chemically Competent Cells are provided with supercoiled pUC19 DNA at a concentration of 10 pg/ μ l—use 1 μ l for transformation.

Preparation for Transformation

E. coli EXPRESS Chemically Competent Cells are provided in aliquots of 80 μ l sufficient for two transformation reactions, respectively, of 40 μ l each.

Transformation is performed by heat shock at 42°C, followed by incubation on ice.

To ensure successful transformation results, the following precautions must be taken:

- For best results, the ligation reaction must be purified or heat killed at 70°C for 15 minutes before transformation.
- All microcentrifuge tubes must be thoroughly pre-chilled on ice before use.
- The cells must be completely thawed **on ice** before use.
- For highest transformation efficiency, use the provided Expression Recovery Medium to resuspend the cells after transformation.

E. cloni[®] EXPRESS Chemically Competent Cells

Transformation Protocol

1. Chill sterile culture tubes on ice (17 mm x 100 mm tubes, one tube for each transformation reaction).
2. Remove *E. cloni* EXPRESS cells from the -80°C freezer and thaw completely on wet ice (10-20 minutes).
3. Add 40 µl of *E. cloni* EXPRESS cells to the chilled culture tube.
4. Add 1 µl of ligation reaction or DNA sample to the 40 µl of cells on ice. (Failure to purify or heat-inactivate the ligation reaction may prevent transformation.) Stir briefly with pipet tip; **do not** pipet up and down to mix, which can introduce air bubbles and warm the cells.
5. Incubate on ice for 30 minutes.
6. Heat shock cells by placing them in a 42°C water bath for 45 seconds.
7. Return the cells to ice for 2 minutes.
8. Add 960 µl of room temperature Expression Recovery Medium to the cells in the culture tube.
9. Place the tubes in a shaking incubator at 250 rpm for 1 hour at 37°C.
10. Plate up to 100 µl of transformed cells on LB agar plates containing the appropriate antibiotic.
11. Incubate the plates overnight at 37°C.
12. Transformed clones can be further grown in LB or any other lactose minus medium.

Sample Induction Protocol

1. Inoculate a single colony from a freshly streaked plate into 5 ml of LB medium containing the appropriate antibiotic for the plasmid and host strain.
2. Incubate with shaking at 37°C overnight. To minimize the expression of the target protein prior to induction, add glucose to the growth medium to a final concentration of 0.2% (w/v).
3. Inoculate 50 ml of LB medium containing the appropriate antibiotic with 0.5 ml of the overnight culture prepared in step 2.
4. Incubate with shaking at 37°C until the OD₆₀₀ reaches 0.6 - 0.8.
5. Add IPTG to a final concentration of 1 mM (Prepare a 1 M solution of IPTG by dissolving 2.38 g of IPTG in water and adjust the final volume to 10 ml. Filter sterilize before use). To determine the optimal concentration of IPTG for maximum expression of the target protein test a range of IPTG concentrations from 0.25 - 2 mM.
6. Incubate at 37°C for 3-4 hours. The optimal time for induction of the target protein may vary from 2-16 hours.
7. Place the culture on ice for 10 minutes. Harvest cells by centrifugation at 5,000 x g for 10 minutes at 4°C.
8. Remove the supernatant and store the cell pellet at -20°C (storage at lower temperatures is also acceptable).

E. coli[®] EXPRESS Chemically Competent Cells

Media Recipes

LB Culture Medium for Growth of Transformants

Per liter: 5 g yeast extract
 8 g tryptone
 5 g NaCl

Add all components except glycerol to deionized water. Adjust pH to 7.0 with NaOH. Autoclave and cool to 55°C.

LB Agar Plates

Per liter: 5 g yeast extract
 8 g tryptone
 5 g NaCl
 15 g agar

Add deionized water to 1 liter. Adjust pH to 7.0 with NaOH. Autoclave. Cool to 55°C and add the appropriate filter-sterilized antibiotic (e.g., 30-50 mg kanamycin for kanamycin-resistant transformants; 100 mg ampicillin or carbenicillin for ampicillin-resistant transformants).

For blue/white screening, add 3 ml 100mM IPTG and 10 ml 2% X-gal to the molten agar at 55°C before pouring.

Pour approximately 25 ml per petri plate.

Related Lucigen Products

- Expresso™ T7 Cloning & Expression System
- OverExpress™ Competent Cells
- CloneSmart® Blunt Cloning Kit
- DNATerminator® End Repair Kit
- PCRTerminator® End Repair Kit
- UltraClone™ DNA Ligation & Transformation Kit
- CloneDirect™ Rapid Ligation Kit
- PCR-SMART™ Cloning Kit
- ClonePlex® Library Construction Kit
- pEZSeq™ Blunt Cloning Kit
- cSMART™ cDNA Cloning Kit
- *E. coli*[®] 10G Competent Cells