

# 10 mM dNTP Set, PCR Grade

## Technical Specifications



Catalog No.: 30029-1 750 µl of each of dATP, dTTP, dGTP and dCTP  
30029-2 5 X Cat. No. 30029-1

### Store at -20°C.

For *In Vitro* Research Use Only.  
Not for Drug or Diagnostic use. Not for use in humans or animals.

Lucigen Corporation  
2120 W. Greenview Dr.  
Middleton, WI 53562 USA  
Tel. 608-831-9011  
Fax. 608-831-9012  
Web [www.lucigen.com](http://www.lucigen.com)  
Email [lucigen@lucigen.com](mailto:lucigen@lucigen.com)

<b>Product Description</b>	10 mM dNTP Set, PCR Grade, contains separate vials of 10 mM sodium salt solutions of each of dATP, dCTP, dGTP and dTTP (pH 8.3).
<b>Storage and usage conditions</b>	Store the product at -20°C. Avoid frequent freeze and thaw cycles. Mix well prior to use.
<b>Stability</b>	10 mM dNTP Set, PCR Grade, is stable for one year from the date received if stored at -20°C.
<b>Recommended Reaction Conditions</b>	200 µM each dNTP; 1X Reaction Buffer; 1 µM primers and 1 - 2.5 U Thermostable DNA Polymerase.
<b>Applications</b>	PCR, RT-PCR, Reverse Transcription, DNA labeling reactions and Sequencing/Cycle Sequencing.
<b>Absence of Endonuclease or Nicking Activity</b>	Incubation of 20 µl of each 10 mM dNTP, PCR Grade, with 1 µg of supercoiled pBR322 DNA for 16 hours at 37°C resulted in no detectable conversion to relaxed or linear forms by agarose gel electrophoresis.
<b>Absence of Exonuclease Activity</b>	Incubation of 20 µl of each 10 mM dNTP, PCR Grade, with 1 µg of HindIII-cut lambda DNA for 16 hours at 37°C resulted in no smearing of bands on agarose gels.
<b>Absence of Ribonuclease Activity</b>	Incubation of 20 µl of each 10 mM dNTP, PCR Grade, with fluorescent labeled RNA substrate resulted in no detectable RNase activity.
<b>Quality Control</b>	The 10 mM dNTP Set, PCR Grade, is tested in DNA amplification using a variety of templates and primers.
<b>Purity</b>	>99% pure.

### PLEASE NOTE

Some applications in which Lucigen's 10 mM dNTP Set, PCR Grade, can be used may be covered by patents issued and applicable in the United States and other countries. Because purchase of this product does not include a license to perform any patented application, users of this product may be required to obtain a patent license depending upon the particular application in which the product is used. The PCR process is the subject of European Patent Nos. 201,184 and 200,262 owned by Hoffman-LaRoche. Those patents expired on March 28, 2006. The corresponding PCR process patents in the United States expired on March 29, 2005. It is the sole responsibility of the buyer to ensure that use of the product does not infringe the patent rights of third parties.

### Warranty

Lucigen warrants that its products will conform to the standards stated in its product specification sheets in effect at the time of shipment. Lucigen will replace, free of charge, any product that does not conform to the specifications. This warranty limits Lucigen's liability only to the replacement of the product. THIS WARRANTY IS EXCLUSIVE, AND LUCIGEN MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The stated express warranties, and the remedy provided for breach thereof, are in lieu of all other liability or obligations of Lucigen for any damages whatsoever arising out of or in connection with the delivery, use, performance, or the inability to use any of its products. IN NO EVENT SHALL LUCIGEN BE LIABLE UNDER ANY LEGAL THEORY (INCLUDING BUT NOT LIMITED TO CONTRACT, NEGLIGENCE, STRICT LIABILITY IN TORT, OR WARRANTY OF ANY KIND) FOR ANY INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES (INCLUDING BUT NOT LIMITED TO LOST PROFITS). Without limiting the effect of the preceding sentence, Lucigen's maximum liability, if any, shall not exceed the purchase price paid by buyer for the product.

# 10 mM dNTP Set, PCR Grade

## PCR Setup Protocol

1) PCR amplification is performed by adding Template DNA (10-50 ng of plasmid DNA; 50-200 ng of genomic DNA), reaction buffer, Forward Primer (100 pmol/ $\mu$ l) and Reverse Primer (100 pmol/ $\mu$ l), dNTPs, thermostable DNA polymerase and water.

2) Reaction Setup. Set up PCR amplifications of the desired size (on ice for best results), according to the following chart:

	25 $\mu$ l Reaction	50 $\mu$ l Reaction	100 $\mu$ l Reaction	Final concentration
DNA template (10 ng/ $\mu$ l)	1.0 $\mu$ l	1.0 $\mu$ l	1.0 $\mu$ l	< 50 ng
10 X Reaction Buffer containing 15 mM MgCl <sub>2</sub>	2.5 $\mu$ l	5.0 $\mu$ l	10.0 $\mu$ l	1 X
Forward Primer (100 pmol/ $\mu$ l)	0.25 $\mu$ l	0.5 $\mu$ l	1.0 $\mu$ l	1 pmol/ $\mu$ l
Reverse Primer (100 pmol/ $\mu$ l)	0.25 $\mu$ l	0.5 $\mu$ l	1.0 $\mu$ l	1 pmol/ $\mu$ l
10 mM dNTP Set (add indicated amount for each dNTP)	0.5 $\mu$ l	1.0 $\mu$ l	2.0 $\mu$ l	0.2 mM each
Thermostable DNA Polymerase (5U/ $\mu$ l)	0.5 $\mu$ l	0.5 $\mu$ l	0.5 $\mu$ l	2.5U
Water, Nuclease-free	18.5 $\mu$ l	38.5 $\mu$ l	78.5 $\mu$ l	---

3) Gently mix the PCR components in a thin-walled reaction tube and spin briefly in a microcentrifuge. Add a drop of mineral oil if the thermocycler does not have a heated lid.

### 4) PCR Cycling Conditions

**Pre-heat the thermocycler to 94°C.**

Cycling step	Temperature	Time	# of Cycles
Initial Denaturation	94°C	2 min	1
Denaturation	94°C	15-30 sec	25 – 35
Annealing*	50-65°C	15-30 sec	
Extension	72°C	1 min/kb	
Final Extension	72°C	5 -10 min	1
Hold	4°C	Indefinitely	1

**\*Anneal at T<sub>m</sub> of primer  $\pm$  2°C.**

5) After completion of the PCR, a 5- $\mu$ l aliquot of the reaction is loaded onto an agarose gel for analysis or size selection.