

# ***Endura***<sup>™</sup> **Chemically Competent and Electrocompetent Cells**

**IMPORTANT!  
-80°C Storage Required  
Immediately Upon Receipt**

**Lucigen<sup>®</sup> Corporation**  
Advanced Products for Molecular Biology

2120 W. Greenview Drive  
Middleton, WI 53562 USA  
Toll Free: (888) 575-9695  
Phone: (608) 831-9011  
FAX: (608) 831-9012  
E-mail: [lucigen@lucigen.com](mailto:lucigen@lucigen.com)  
Internet: [www.lucigen.com](http://www.lucigen.com)

# Endura™ Competent Cells

## Contents

Components & Storage Conditions.....	2
Endura Electrocompetent Cells.....	3
Endura Electrocompetent Cells Transformation Protocol.....	3
Endura Chemically Competent Cells .....	4
Endura Chemically Competent Cells Transformation Protocol .....	5
Related Products.....	6

## Components & Storage Conditions

Lucigen's EnduraCompetent Cells are packed in Lucigen's DUO format for your convenience:

### Electrocompetent cells:

**DUO:** 50 µL per vial, enough for two 25-µL reactions

### Chemically competent cells:

**DUO:** 80 µL per vial, enough for two 40-µL reactions

**SOLO:** 40 µL per vial, enough for one reaction

The cells are shipped on dry ice in one container, along with Recovery Medium and supercoiled control pUC19 DNA at 10 pg/µl. Please refer to Table 1 for a complete listing of *Endura* Electrocompetent Cells, efficiencies, and catalog numbers.

All Competent Cells require storage at -80° C.

**Table 1: *Endura* Competent Cells available from Lucigen**

STRAIN	Format	Efficiency (cfu/µg pUC19)	Transformations	Catalog #	Storage
<i>Endura</i> Chemically Competent Cells	SOLO	> 1 x 10 <sup>8</sup>	12 ( 12 x 40 µL ) 24 ( 24 x 40 µL )	60241-1 60241-2	-80°C
	DUO	> 1 x 10 <sup>8</sup>	4 ( 2 x 80 µL ) 12 ( 6 x 80 µL ) 24 ( 12 x 80 µL )	60240-0 60240-1 60240-2	-80°C
<i>Endura</i> Electrocompetent Cells	DUO	> 1 x 10 <sup>10</sup>	4 ( 2 x 50 µL ) 12 ( 6 x 50 µL ) 24 ( 12 x 50 µL )	60242-0 60242-1 60242-2	-80°C
Recovery Medium	-		4 ( 1 x 12 mL ) 12 ( 1 x 12 mL ) 24 ( 2 x 12 mL ) ( 8 x 12 mL )	---- ---- ---- 80026-1	-80 to +20°C
*Supercoiled pUC19 DNA (10 pg/µL)	-		( 1 x 20 µL )	----	-80 to -20°C

\*Provided as a positive control for transformation.

# Endura™ Competent Cells

## Endura Competent Cells

Lucigen's Endura Competent Cells are noted for their ability to clone and maintain plasmid DNA sequences that are otherwise unstable, such as those with inverted repeats or sequences that are otherwise prone to undesired recombination events.

Such unstable inserts are commonly found in retroviral DNA, and can prove exceptionally difficult to isolate and propagate in traditional cloning strains. Lucigen's Endura Competent Cells are available in high efficiency chemically competent or electrocompetent versions to make difficult cloning work more reliable and efficient.

### Transformation Efficiency:

Chemically Competent Cells:  $\geq 1 \times 10^8$  cfu/ $\mu$ g  
Electrocompetent Cells:  $\geq 1 \times 10^{10}$  cfu/ $\mu$ g

### Genotype

*recA13 supE44 ara-14 galK2 lacY1 proA2 rpsL20(Str<sup>R</sup>) xyl-5 λ<sup>-</sup> leu mtl-1 F<sup>-</sup> mcrB mrr  
hsdS20(r<sub>B</sub><sup>-</sup>, m<sub>B</sub><sup>-</sup>)*

## Usage Guidelines: Electrocompetent Cells

### Preparation for Transformation

Transformation is carried out in a 0.1 cm gap cuvette using 25  $\mu$ L of *Endura* Electrocompetent Cells. Optimal settings for electroporation are listed in the table below. Typical time constants are 3.5 to 4.5 msec.

Optimal Setting	Alternate Settings (~ 20-50% lower efficiencies)
1.0 mm cuvette 10 $\mu$ F 600 Ohms 1800 Volts	1.0 mm cuvette 25 $\mu$ F 200 Ohms 1400 – 1600 Volts

### Suggested Electroporation Systems:

Bio-Rad Micro Pulser #165-2100; Bio-Rad E. coli Pulser #165-2102; Bio-Rad Gene Pulser II #165-2105; BTX ECM630 Electroporation System; Eppendorf Electroporator 2510.

Optimal transformation control reactions include electroporation with 1  $\mu$ L (10 pg) of supercoiled pUC19 DNA.

To ensure successful transformation results, the following precautions must be taken:

- Ligation reactions performed with Lucigen's CloneDirect™ Ligation Buffer (included with Lucigen's Cloning or Ligation Kits) must be heat killed at 70°C for 15 minutes before transformation. The ligation reaction can be used directly after heat inactivation, without purification of the ligation products.
- DNA samples in other buffers must be purified and dissolved in water or a buffer with low ionic strength, such as TE. The presence of salt in the electroporation sample leads to arcing at high voltage, resulting in loss of the cells and DNA.

# Endura™ Competent Cells

---

- Microcentrifuge tubes and electroporation cuvettes must be thoroughly pre-chilled on ice before use. Successful results can be obtained with cuvettes from many sources including BTX (Model 610), Eppendorf (Cat. #940001005), and BioRad (Cat. #165-2089).
- The cells must be completely thawed **on ice** before use.
- For highest transformation efficiency, use the provided Recovery Medium to resuspend the cells after electroporation. Use of SOC or other media will result in lower transformation efficiencies.
- Prepare nutrient agar plus appropriate antibiotic.

## Transformation Protocol

Use the following protocol for cells provided in microfuge tubes.

1. Have Recovery Medium and 17 mm x 100 mm sterile culture tubes readily available at room temperature (one tube for each transformation reaction). Transformation efficiency may decrease with the use of SOC or other media.
2. Place electroporation cuvettes (0.1 cm gap) and microcentrifuge tubes on ice (one cuvette and one microfuge tube for each transformation reaction).
3. Remove cells from the -80°C freezer and place on wet ice until they thaw **completely** (10-20 minutes).
4. When cells are thawed, mix them by tapping gently. Aliquot 25 µL of cells to the chilled microcentrifuge tubes on ice. (Omit this step if using *E. coli* SOLOs, which contain 25 µL of cells per tube).
5. If using ligation buffer from any Lucigen Cloning Kit, add 1 µL of the heat-denatured ligation reaction to the 25 µL of cells on ice. Failure to heat-inactivate the ligation reaction will prevent transformation. Stir briefly with pipet tip; **do not** pipet up and down to mix, which can introduce air bubbles and warm the cells. Use of more than 2 µL of ligation mix may cause electrical arcing during electroporation.

For ligation reactions using other commercial kits, please refer to the manufacturer's instructions.

6. Carefully pipet 25 µL of the cell/DNA mixture into a chilled electroporation cuvette without introducing bubbles. Quickly flick the cuvette downward with your wrist to deposit the cells across the bottom of the well. Electroporate according to the conditions recommended above.
7. Within 10 seconds of the pulse, add 975 µL of Recovery Medium to the cuvette and pipet up and down three times to resuspend the cells. Transfer the cells and Recovery Medium to a culture tube.
8. Place the tube in a shaking incubator at 250 rpm for 1 hour at 37°C.
9. Spread up to 100 µl of transformed cells on nutrient agar plates containing the appropriate antibiotic.
10. Incubate the plates overnight at 37°C.
11. Transformed clones can be further grown in any rich culture medium.

## Usage Guidelines: Chemically Competent Cells

### Preparation for Transformation

To ensure successful transformation results, the following precautions must be taken:

- For best results, Lucigen CloneSmart® ligation reactions must be heat killed at 70°C for 15 minutes before transformation. Alternately, the reactions may be purified, if desired. For other ligation reactions, follow the manufacturer's recommendations.
- Prepare nutrient agar plus antibiotic for selection.
- All microcentrifuge tubes must be thoroughly pre-chilled on ice before use.
- The cells must be completely thawed **on ice** before use.
- For highest transformation efficiency, use the provided Recovery Medium to resuspend the cells after transformation.

### Transformation Protocol

1. Prepare nutrient agar (e.g., LB) plates with antibiotic for selection.
2. Chill sterile culture tubes on ice (17 mm x 100 mm tubes, one tube for each transformation reaction).
3. Remove cells from the -80°C freezer and thaw completely on wet ice (10-20 minutes).
4. Add 40 µL of cells to the chilled culture tube.
5. Add 1-4 µL of heat-inactivated ligation reaction or DNA sample to the 40 µL of cells on ice. (Failure to heat-inactivate—70°C for 15 minutes—or otherwise purify, the ligation reaction may prevent transformation.) Stir briefly with pipet tip; **do not** pipet up and down to mix, which can introduce air bubbles and warm the cells.
6. Incubate on ice for 30 minutes.
7. Heat shock cells by placing them in a 42°C water bath for 45 seconds.
8. Return the cells to ice for 2 minutes.
9. Add 960 µL of room temperature Recovery Medium to the cells in the culture tube. When using these cells with a cloning kit, follow the Recovery Medium volume given in that kit manual.
10. Place the tubes in a shaking incubator at 250 rpm for 1 hour at 37 °C.
11. Plate up to 100 µL of transformed cells on nutrient agar plates containing the appropriate antibiotic.  
Note: the quality of LB plates varies widely. Transformants plated on LB may grow slowly.
12. Incubate the plates overnight at 37°C.
13. Transformed clones can be further grown in any rich culture medium (e.g., LB, or TB).

## Related Lucigen Products

- BigEasy® Linear Cloning Kit
- CloneSmart® Blunt Cloning Kits
- DNATerminator® End Repair Kit
- PCRTerminator® End Repair Kit
- UltraClone™ DNA Ligation & Transformation Kits
- CloneDirect™ Rapid Ligation Kit
- PCR-SMART™ Cloning Kits
- ClonePlex® Library Construction Kit
- pEZSeq™ Blunt Cloning Kits
- cSMART™ cDNA Cloning Kits
- *E. coli*® EXPRESS Electrocompetent Cells
- OverExpress™ Competent Cells

## Notice of Limited Label License, Copyright, Patents, Warranties, Disclaimers, and Trademarks

Copyright© 2001-2010 by Lucigen Corp. All rights reserved. Lucigen, BigEasy, ClonePlex, CloneSmart, CopyRight, DNATerminator, *E. coli*, EconoTaq, ExCyto, Lucigen, NanoClone, PCR Terminator, pJAZZ, pSMART, pTrueBlue, and PyroPhage are registered trademarks of Lucigen Corporation. bSMART, Blue/White cloning without the blues, ChimeraFree, CloneDirect, cSMART, eLucidations, GapFree, PCR SMART, pEZ, pEZSeq, pGC, pLEXX-AK, pRANGER, Replicator, Softag, The Cloning Company, The Molecular Cloning Company, Transformance, TSA, and UltraClone are trademarks of Lucigen Corporation. Hydroshear and GeneMachines are trademarks of Genomic Solutions, Inc. DH10B is a trademark of Invitrogen Corp.

Lucigen's products are sold for research use only and are not to be used in humans or for medical diagnostics. Lucigen's liability with respect to any of its products is limited to the replacement of the product. No other warranties of any kind, expressed or implied, including without limitation, any implied fitness for any particular use, are provided by Lucigen. Lucigen is not liable for any direct, indirect, incidental or consequential damages arising out of or in connection with the use or inability to use any of its products.

If the purchaser is not willing to accept these use limitations, Lucigen Corporation is willing to accept return of the product for a full refund. For information on obtaining a license, contact Lucigen Corporation, 2120 W. Greenview Dr., Middleton, WI 53562. Email: [Lucigen@lucigen.com](mailto:Lucigen@lucigen.com). Phone: 608-831-9011. Fax 608-831-9012.